

Espec Software User Guide

Introduction

Welcome to Espec, a pioneering software designed to enhance the functionality of second-generation spectrophotometer hardware. Crafted with the core principles of the iGEM spirit, Espec aims to streamline and simplify spectrophotometry for both novices and experts.

Prerequisites

Software License: Espec is open-source software, released under the MIT License. Copyright (c) 2023 by CSMU-Taiwan.

Required Hardware: Second-generation spectrophotometer setup as detailed in the corresponding hardware assembly guide.

Getting Started

1. Initiating Espec:

- Launch the Espec application.

2. Connection to Arduino:

- Enter the IP address of the ESP32.
- The software will validate the IP and provide options to reconnect if any discrepancies arise.

3. LED Customization:

- Users can select from red, green, or blue LEDs based on their experimental needs.

4. Data Harvesting from Arduino:

- The software retrieves data every 5 seconds.
- Data categories displayed include illumination, blank, test, and computational values.


5. Calibration Curve Insights:

- Users can add, save, or delete data points.
- Options to clear all data or generate a trendline for in-depth analysis are available.

6. Trendline Data Exposition:

- The software displays predicted concentrations, trendline equations, and R-squared values for comprehensive analysis.

Espec Software Operational Flowchart



Navigating the User Interface

1. **Dashboard:** This central hub allows users to connect to Arduino, choose LED colors, and view real-time data.

Spectrometer

Back

Connect to Arduino

192.168.0.224
Confirm
Reconnect

Adjust LED Color

Red LED
Green LED
Blue LED

Data From Arduino

Data Type	Value	Actions
Lux Value	319	
Blank Value	317	Blank
Test Value	321	Test
Calculate Value	1.0126	

2. Data Visualization and Trendline

Analysis: Transform raw data into insightful graphs and trendlines.

Customize, save, or reset visual representations based on user preferences.

Limitations & Future Prospects

While Espec offers a plethora of features, it's primarily optimized for second-generation hardware. Future updates seek to expand device compatibilities and integrate more advanced features based on user feedback.

Conclusion & Acknowledgments

Espec is a testament to the potential of open-source solutions in synthetic biology research. As we embark on this journey, we extend our gratitude to the iGEM community for their continuous support and inspiration. We invite everyone to test Espec and contribute to its evolution.

Acknowledgments:

This software was developed by the 2023 CSMU-Taiwan iGEM team. We extend our gratitude to our team members, advisors, and the broader iGEM community for their invaluable feedback and support.

